



Imagerie de Diffusion IVIM

Alain Luciani, Frederic Pigneur, Emmanuel Itti, Alain Rahmouni

CHU Henri Mondor, Université Paris Est Creteil

alain.luciani@aphp.fr

Comment je fais ?

Antennes
Séquences

IVIM DWI-MRI

- Pré-requis technologiques et instrumentaux





Applications Cliniques

Comment je fais ?


Antennes

Surface

Antenne Corps



- Comment gérer les mouvements ?
 - TE minimum
 - Imagerie parallèle
 - Asservissement respiratoire ?



12 volontaires sains

Trois méthodes de gating

- Respiration libre
- Gating Respiratoire
- Gating Cardiaque

Etude DWI IVIM : 9 valeurs de b : 0, 30, 60, 100, 150, 200, 400, 600, et 900 sec/mm²

T

| Parameter | FB DWI | ET DWI | RT DWI |
|---------------------------------|-----------|-----------------|-----------------|
| Echo time (msec) | 60 | 60 | 60 |
| Repetition time (msec) | 2100 | 3243–7142* | 2100 |
| Echo-planar imaging factor | 115 | 115 | 115 |
| Receiver bandwidth (Hz) | 1594 | 1594 | 1594 |
| Field of view (mm) | 340 × 256 | 340 × 256 | 340 × 256 |
| Matrix size | 196 × 115 | 196 × 115 | 196 × 115 |
| No. of signals acquired | 4 | 4 | 4 |
| Section thickness (mm)/gap (mm) | 7/1.4 | 7/1.4 | 7/1.4 |
| No. of sections acquired | 20 | 20 | 20 |
| Acquisition time (min) | 3.7 | 8.3 (6.5–13.0)† | 6.8 (6.4–10.8)† |


12 volontaires sains


Trois méthodes de gating

- Respiration libre
- Gating Respiratoire
- Gating Cardiaque

Etude DWI IVIM : 9 valeurs de b : 0, 30, 60, 100, 150, 200, 400, 600, et 900 sec/mm^2

T





| Parameters | Within-Subject CV (%) | | | P-values* |
|------------|------------------------|-----------------------|------------------------|------------------|
| | FB DWI | ET DWI | FI DWI | |
| ADC | | | | |
| Right lobe | 4.33 (2.95, 5.70) | 3.17 (1.79, 4.55) | 2.33 (1.50, 3.16) | .279, .022, .334 |
| Left lobe | 4.74 (3.12, 6.37) | 4.12 (2.37, 5.88) | 6.96 (2.16, 11.75) | .658, .365, .245 |
| Mean | 4.15 (2.66, 5.64) | 3.40 (1.75, 5.05) | 4.59 (1.24, 7.93) | .568, .806, .546 |
| D-value | | | | |
| Right lobe | 4.11 (2.51, 5.71) | 4.05 (2.81, 5.30) | 3.19 (1.83, 4.54) | .960, .378, .282 |
| Left lobe | 12.51 (4.61, 20.40) | 5.34 (3.46, 7.21) | 16.17 (7.35, 25.00) | .087, .072, .023 |
| Mean | 7.57 (3.45, 11.68) | 4.52 (2.96, 6.08) | 9.73 (4.53, 14.93) | .194, .097, .074 |
| F-value | | | | |
| Right lobe | 15.77 (10.63, 20.90) | 9.86 (7.13, 12.59) | 11.85 (5.91, 17.78) | .004, .306, .525 |
| Left lobe | 16.81 (9.60, 24.02) | 8.47 (6.79, 10.74) | 12.10 (5.22, 18.98) | .002, .003, .347 |
| Mean | 13.83 (9.49, 18.18) | 7.60 (5.79, 9.42) | 10.08 (4.91, 15.19) | .002, .151, .348 |
| D* value | | | | |
| Right lobe | 74.16 (43.69, 104.63) | 71.87 (47.53, 96.21) | 75.03 (55.04, 95.03) | .919, .964, .813 |
| Left lobe | 156.61 (96.72, 216.49) | 88.34 (31.57, 145.11) | 103.75 (49.18, 158.33) | .119, .224, .181 |
| Mean | 95.60 (56.31, 134.90) | 57.05 (40.00, 74.10) | 63.74 (27.90, 99.58) | .103, .208, .756 |

- Respiration : Impact sur WB – imaging ?
 - 10 volontaires
 - DWI SE EPI sequence
 - TE 76, b0 , 400s/mm2
 - Respiration libre

Short- and Midterm Reproducibility of Apparent Diffusion Coefficient Measurements at 3.0-T Diffusion-weighted Imaging of the Abdomen¹

Purpose: To test the hypothesis that there is no significant variability in apparent diffusion coefficients (ADCs) at assessment of the short- and midterm reproducibility of ADC measurements in a healthy population.

Materials and Twenty healthy male volunteers were enrolled in this pro-

Adam C. Brailwaile, MD
Brian M. Dale, PhD
Daniel T. Boll, MD
Elmar M. Merkle, MD

SIGNAL RESEARCH ■ GASTROINTESTINAL IMAGING





Figure 1: (a) Axial two-dimensional DW image of liver, spleen, and pancreatic tail, acquired with b value of 0 sec/mm². Three ROIs (1–9) each are placed in these anatomic areas. (b) The same ROIs have been copied from the DW image onto the corresponding ADC map.

- Variation intra-individuelle ADC $\pm 14\%$
- Seuil de signification $> 27\%$



- Respiration : Impact sur WB-DWI ?

JOURNAL OF MAGNETIC RESONANCE IMAGING 30:561-568 (2009)

Original Research

**Diffusion-Weighted Imaging of the Liver:
Comparison of Navigator Triggered and Breathhold
Acquisitions**

Bachir Taouli, MD,^{1*} Alison Sandberg, MD,¹ Alto Stemmer, PhD,² Tejas Parikh, MD,¹ Samson Wong, MD,¹ Jian Xu, PhD,³ and Vivian S. Lee, MD, PhD¹






Figure 2. A 75-year-old man with hemangioma of the right posterior hepatic lobe (long arrows) demonstrated with breathhold DWI and respiratory triggered (PACE) DWI for b = 0, 50, and 500 s/mm² and corresponding ADC maps. An associated left renal cyst is present (short arrow). PACE DWI shows better image quality, and better lesion conspicuity than BH DWI. ADC map for PACE DWI shows better lesion delineation and less noise contamination.


- Nombre de gradients de diffusion ?

Mêmes paramètres ; Nex=2







b0




b10




b20




b30




b40




b100



b200





b400



b800

- Choix des valeurs de b ?



- Choix des facteurs de diffusion b – **Détection**





Figure 4. Mean (\pm SD) metastasis to liver CNR of 35 metastases larger than 10 mm detected with IVIM-DWI. The maximum CNR was obtained with $b = 100 \text{ s/mm}^2$.

- Principes de la diffusion



Diffusion moléculaire



Composante perfusive

• Choix des valeurs de b ?



- b 0, 500, 1000 s/mm²
- Modèle animal CCl4

JOURNAL OF MAGNETIC RESONANCE IMAGING 25:122-128 (2007)



Original Research

Assessment of Diffusion-Weighted MR Imaging in Liver Fibrosis


Laurence Annet, MD, PhD,^{1*} Frank Peeters, PhD,¹ Jorge Abarca-Quinones, BS,¹ Isabelle Leclercq, MD, PhD,² Pierre Moulin, MD,³ and Bernard E. Van Beers, MD, PhD¹



- Choix du facteur b - **Quantification**



- Choix du facteur b - **Quantification**



$$S = S_0[(1-f) \cdot \exp(-b \cdot D) + f \cdot \exp(-b(D+D^*))]$$

- Principes de la diffusion



$$B = \frac{S(b)}{S(0)} = \exp(-bD).F = \underbrace{f.\exp(-b(D + D^*))}_{\text{Microcirculation Composante perfusive}} + \underbrace{(1 - f).\exp(-bD)}_{\text{Diffusion moléculaire}}$$


$b < 100 \text{ s/mm}^2$: composante perfusive non négligeable

$b > 100 \text{ s/mm}^2$: composante perfusive négligeable

Microcirculation
Composante
perfusive

Diffusion
moléculaire

• Application clinique



Alain Luciani, MD, PhD
Alexandre Vignaud, PhD
Madeleine Casat, MD
Jeanne Truc Van Nhieu, MD
Ariane Multat, MD, PhD
Lucile Ruel, PhD
Alexis Laurent, MD, PhD
Jean-François Daux, MD, PhD
Pierre Bruguier, MD
Alain Rahmouni, MD

Liver Cirrhosis: Intravoxel Incoherent Motion MR Imaging—Pilot Study¹

Purpose: To retrospectively evaluate a respiratory-triggered diffusion-weighted (DW) magnetic resonance (MR) imaging sequence combined with parallel acquisition to allow the calculation of pure molecular-based (D) and perfusion-related (D^* , f) diffusion parameters, on the basis of the intravoxel incoherent motion (IVIM) theory, to determine if these parameters differ between patients with cirrhosis and patients without liver fibrosis.

Materials and Methods: The institutional review board approved this retrospective study; informed consent was waived. IVIM DW imaging

ORIGINAL RESEARCH ■ GASTROINTESTINAL IMAGING

Table 4

Mean global ADC and D , D^* , and f Parameters Assessed in 37 Patients in Pilot Study

| Parameter | Healthy Liver Group | Cirrhotic Liver Group | P Value |
|--|---------------------|-----------------------|----------|
| ADC ($\times 10^{-3}$ mm ² /sec) | 1.39 \pm 0.2 | 1.23 \pm 0.4 | .03 |
| D ($\times 10^{-3}$ mm ² /sec) | 1.10 \pm 0.7 | 1.19 \pm 0.5 | .4 (NS) |
| D^* ($\times 10^{-3}$ mm ² /sec) | 79.1 \pm 18.1 | 59.4 \pm 20.0 | .008 |
| f (%) | 27.0 \pm 5.3 | 30.0 \pm 5.7 | .07 (NS) |

Note.—Unless otherwise indicated, data are means \pm standard deviations. NS = not significant.


IVIM €

- Modèle animal murin
- TAA 0.03%
- Eau de boisson
- Adaptation poids 6 semaines




Saci
anat


ADC / Fibrose




D* / Fibrose




/BK/6J




- Choix du facteur b



- T2 WI
- Limite effet « perfusion »
- Filtre vitesse



$b = 0$ s/mm²



$b = 50$ s/mm²

Pourquoi améliorer analyse DWI ?

RÉPARTITION DES TUMEURS HÉPATIQUES

| Tumeur N = 140 | Nombre | | Taille (mm) | |
|-------------------|--------|------|-------------|---------------------|
| | Total | Naïf | Moyenne | Déviati on standard |
| CHC | 50 | 27 | 52,78 | 44,9 |
| Métastase * | 29 | 7 | 33,69 | 17,22 |
| HNF | 14 | | 49,93 | 25,38 |
| Adénome | 3 | | 34,7 | 28,9 |
| Angiome | 34 | | 30,65 | 16,98 |
| Kyste | 10 | | 32,3 | 15,43 |

Pourquoi améliorer analyse DWI ?

Exemple des Tumeurs Hépatiques






Figure 1 : Paramètres IVIM des lésions hépatiques. ★★★ p < 0,0005 - ** 0,0005 < p < 0,005 - * 0,005 < p < 0,05





Pourquoi améliorer analyse DWI ?

Exemple des Tumeurs Hépatiques







Pourquoi améliorer analyse DWI ?

Exemple des Tumeurs Hépatiques



Pourquoi améliorer analyse DWI ?

Exemple des Tumeurs Hépatiques




Pourquoi améliorer analyse DWI ?

Exemple des Tumeurs pancréatiques

Diffusion weighted images were acquired using a single-shot echo-planar imaging (SE-EPI) pulse sequence in expiratory breath-hold with the following imaging parameters: TR=1300 ms, TE=60 ms, FOV=350×273 mm², matrix size=100×78, 14 slices, slice thickness/gap=5/0.25 mm, spectral fat saturation, 4 averages, bandwidth=3000 Hz/pixel, k-space based parallel imaging technique (GRAPPA) acceleration factor of two, *b*-values=0, 50, 100, 150, 200, 300, 400, 600, and 800 s/mm² and a total measurement time of 12 minutes. The

Table 1 Comparison of mean IVIM-parameters and ADC in AIP-examinations

| | AIP initial | AIP 1 st follow up | AIP 2 nd follow up | p-value |
|--|-------------|-------------------------------|-------------------------------|----------|
| mean <i>f</i> (%) | 10.5±4.3 | 17.1±6.7 | 20.4±4.3 | p<0.0001 |
| mean <i>D</i> (*10 ⁻³ mm ² /s) | 1.19±0.28 | 1.13±0.22 | 1.13±0.2 | p=0.668 |
| mean <i>D</i> * (*10 ⁻³ mm ² /s) | 52.9±59.3 | 34.1±51.1 | 14.0±4.4 | p=0.412 |
| mean ADC (*10 ⁻³ mm ² /s) | 1.28±0.18 | 1.34±0.25 | 1.38±0.24 | p=0.549 |



DOI: 10.1007/s00130-015-0961-4

HEPATOBIILIARY-PANCREAS

IVIM DW-MRI of autoimmune pancreatitis: therapy monitoring and differentiation from pancreatic cancer

Miriam Klauß¹ · Klaus Maier-Hein² · Christine Tjaden³ · Thilo Hackert⁴ · Lars Grenacher¹ · Bram Stieltjes⁴

- Nouvelles approches
 - Etudes non gaussiennes : Kurtosis
 - Membranes / Barrières




Fig. 1. Plot of the logarithm of signal decay vs. b value in a selected WM ROI belonging to the corpus callosum, averaged over all pixels. The dashed line represents the predicted mono-exponential decay at $b=1000 \text{ s/mm}^2$, while the dotted line represents the noise floor, calculated as the average of the pixel intensities in a ROI outside the brain.





Figure 2. Comparison of DTI and DKI fitting models. For DTI, the logarithm of diffusion-weighted signal intensity (circles) as a function of the b-value is fit, for small b-values, to a straight line. In brain, this fit is often based on the signal for $b=0$ and $b=1000 \text{ s/mm}^2$. For DKI, the logarithm of the signal intensity is fit, for small b-values, to a parabola. In brain, this fit may be based on the signal for $b=0$, $b=1000$, and $b=2000 \text{ s/mm}^2$.

Recommandations

1. Chronogramme séquence DWI

1. Vérifier suppression de graisse
2. Vérifier monitoring respiratoire
 1. TR et fréquence respiratoire
 2. PACE trigger



Courtesy: Dr Alexandre Vignaud, PhD
UNIRS, I2BM, CEA, Neurospin

Simplified diagram of hepatic Single-Shot Echo-Planar-Imaging sequence (SS EPI).




- Respiratory trigger using PACE (Prospective Acquisition CorrEction) for the respiration.
- During the expiration plateau as many as possible slices are acquired for a given b-value (if it is not enough concatenations are used).
- Fat Saturation is done using a SPAIR (SPectral AdiAbatic Inversion Recovery) preparation and Inversion Time (TI) 230ms (at 3T).

Recommandations


2. Comment optimiser les séquences

1. Choix des valeurs de b
2. Optimiser RSB
3. Compromis...

- Au moins 2 valeurs de b
 - b faible → Détection → 50 s/mm²
 - b élevé → Characterization → 800 s/mm²
- Ou IVIM
 - Multiples valeurs de b values < et > 100 s/mm²



$b = 0 \text{ s/mm}^2$



$b = 50 \text{ s/mm}^2$

Recommandations

2. Comment optimiser les séquences

1. Choix des valeurs de b
2. Optimiser RSB
3. Compromis...

Courtesy: Dr Alexandre Vignaud, PhD
UNIRS, I2BM, CEA, Neurospin

☐ Champ magnétique

☐ 3 T > 1.5T


☐ Antennes

☐ Surface

☐ Canaux multiples

☐ Nb excitations variable pour chaque valeur d b

☐ Acquisitions monopolaires → Réduit TE → augmente SNR



Recommandations

- ❑ Conserver résolution spatiale
- ❑ réduction de durée du train EPI → Réduit distorsion → augmente bande passante

- ❑ Augmenter résolution spatiale → Durée train EPI plus longs → Distorsion
- ❑ Augmentation de bande passante → Réduit RSB

2. Comment optimiser les séquences





Temps acquisition



1. Choix des valeurs de b



2. Optimiser RSB

3. Compromis...





Siemens MRI
Verio 3T

- b 800s/mm²
- 32 channel
- iPAT 2
- ADC maps



Siemens MRI
Skyra 3T

- b 800s/mm²
- 60 Channel
- Caipi parallel
- SMS
- ADC Maps




Messages

- Instrumentation
- Instrumentation
- Instrumentation

Applications cliniques

- Pas encore définies
- Standardisation ?
- Sélectionner études cliniques adaptées....



"AREN'T THERE ENOUGH PROBLEMS IN THE WORLD ALREADY?"